ReFa

ReFa CARAT RAY

GUIDE BOOK

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
English

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ENGLISH

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Thank you for purchasing ReFa CARAT RAY.

ReFa represents a sincere desire to see women harness the beauty that flows from within and without them——.

It is a brand that expresses a wish to see you delight in discovering your true radiance with just a minimum of effort.

It also symbolizes an intention to produce a spontaneous smile on your face as you enjoy a sense of fulfillment from simply using our product.

ReFa is here for your everyday pleasure.

Deep kneading motion.
A sensational way to guide your skin
to brilliant beauty.

* Photos of RF-PC2019B are used
in this document.

■ Esthetic Method

This device features a 360° multi-angle structure that takes into account factors such as roller shape and angle of incidence to replicate the kneading techniques of a professional.

This innovative design easily adapts to complex body contours, virtually ensuring smooth and uninterrupted kneading regardless of the device's tilt angle.

■ Microcurrent

ReFa rollers emit a microcurrent by taking in light through a built-in solar-panel. ReFa CARAT RAY harnesses light with its built-in solar panels to generate mild electrical microcurrent to a higher degree than the previous ReFa CARAT devices.

■ Platinum Coating (RF-PC2019B only)

To protect your skin during contact, the surface of the device's rollers are coated with platinum.


■ Waterproof Construction

Conforms with JIS standards (equivalent to IPX7).
So the rollers can be used in the bathtub.

Rolling Method

The strength in which the roller pulls up your skin can be adjusted by changing the angle at which you hold the handle.


<When Pushing the Rollers Forward>

Gradually raise the end of the handle to check the pull of the rollers.

*Lower the handle if you begin to feel pain.

<When Drawing the Rollers Back>

Gradually lower the end of the handle to check the pull of the rollers.

*Raise the handle if you begin to feel pain.

Tips on Rolling

Roll in the direction indicated by the arrows and increase or decrease pressure in accordance with the thickness of the lines.


Guidance Arrows:

The photos provided in this manual are overlaid with the following arrows to show how and which direction to roll over specific parts of the body for optimum results.

-  : Roll only in one direction
-  : Roll in both directions
-  : Use more force as you roll from the narrow to the wider end.

How To Use

- Refer to the instructions in the following pages and use according to your preference.
- As a general rule, roll back and forth over each body part about 5 to 10 times.

Please be sure to read the instruction manual before use

- When using this device for the first time, roll it gently a few times over an area of the body that is not very noticeable such as your upper arm, and make sure that there are no abnormal skin reactions after some time has elapsed.
- Persons with metal allergies, or who have conditions such as skin infections or dermatitis, should not use the product.
- Women in the early unstable stages of pregnancy or immediately after childbirth, or women who are pregnant and not in good health should not use this product.
- Remove any accessories (necklaces, piercings, etc.) during use, as they may get caught in the rollers.
- Tie your hair up or wear a hair band when using the device so your hair does not get caught.
- Do not use along the spine or in other locations where bones protrude.
- Please do not use anywhere other than as specified in this manual.

Decolletage & Neck

For a beautiful décolletage and a supple neckline.

Roll along the indentation beneath the collar bone.

Roll up and down along the line of the neck between the head and the base of the shoulder and neck.

Face & Head

For a supple, youthful face and a softer scalp.

Consciously and steadily roll along the line of the lower jaw.

Roll along your cheekbones from the corners of your mouth to below your ears.

Roll up and down between the corner of your mouth and the side of your nose. While doing this, lift the cheek upward with one roller.

Position one roller above the outer corner of the eye and roll it along the side of the head.

With your eyes closed, slowly roll from the forehead up to the top of the head.

*Make sure your hair does not get caught in the product.

Roll up and down with one roller positioned above the ear.

*Make sure your hair does not get caught in the product.

Arm

Knead the skin for a beautiful upper arm contour.

Sandwich the flesh between thumb and forefinger with the rollers, and move the roller pressed against the back of the hand in a semicircular motion.

*Do not place fingers individually between the rollers. Doing so may result in injury to the joints of a finger.

Roll from the wrist to the elbow with the inside of the lower arm held between the rollers.

Pressing against the center of the fleshy portion of the palm between the wrist and the little finger, roll towards the little finger.

*Do not place fingers individually between the rollers. Doing so may result in injury to the joints of a finger.

Roll your upper arms toward your armpit.

*This is an area where the skin can be easily pinched by the rollers. Always be careful when adjusting the angle and pressure of the rollers against the skin.

Bust

To achieve the ideal bust line, always move the roller towards the center of the chest.


Pressing down below the armpit, roll towards the center of the breast.

Roll up and down, applying pressure upwards towards the breast.

Waist

Tighten your skin for a feminine beauty.

Roll along the body between the lower bust line and the waistline.

*This is an area where the skin can be easily pinched by the rollers.

Always be careful when adjusting the angle and pressure of the rollers against the skin.

In the abdomen area, draw upwards as you roll diagonally.

*This is an area where the skin can be easily pinched by the rollers.

Always be careful when adjusting the angle and pressure of the rollers against the skin.

Legs

Firmly approach those troublesome areas around the knees and thighs.

Roll over the entire sole of the foot.
Finish up by pressing one roller
against the area indicated by the "o".

Roll over the sides and instep of
the foot. This can be especially
comforting after a bath.
*Do not apply excessive force when
rolling along the joints of a finger.
Doing so can result in injury.

Roll up and down around the ankle.
Apply an amount of pressure that is
most comfortable to you.
*Be careful your Achilles tendon does
not get pinched in the product.

Roll along the calf from the ankle to
behind the knee.
*Be careful your Achilles tendon does
not get pinched in the product.

Roll on the inner side of your thigh.
*This is an area where the skin can be easily pinched by the rollers.
Always be careful when adjusting the angle and pressure of the rollers against the skin.

Roll the area along the back of the thighs between the buttocks and the back of the knees.
*This is an area where the skin can be easily pinched by the rollers.
Always be careful when adjusting the angle and pressure of the rollers against the skin.

For bath-time relaxation

Can be used in the bath as it is waterproof.

*After use, shake it off and wipe with a soft cloth before storing it.

*Images of how to use.

Assisted Rolling


Two people are better than one when it comes to effectively rolling those areas that are hard to reach such as your back or shoulders.

We recommend using on your shoulders and back.

For a smooth and firm neck line

Roll gently along the neck and shoulder.

For an elegant back

Roll on the back, which is hard to reach yourself. Roll up and down along the shoulder blades.

*Do not use along the spine or in other locations where bones protrude.

A Note on Assisted Rolling (RF-PC2019B only)

The hand not holding the device must stay in constant contact with the person's skin.

*Microcurrent will not be circulated if there is no physical contact.

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